

Avid

Axcel

K0856 / K0861

Enhancing life through mobility.



Contact your Territory Sales Manager or email us at sales@meritsusa.com

Scan to learn more!

©Merits Health Products, Inc. rev. 1024

Features

- Independent hydraulic, high-performance suspension
- Large 9" casters to navigate obstacles and rough terrain
- Bright LED lights, included
- Tool-free, rear battery access
- Power recline adjustment 85° - 160°
- Power tilt adjustment 0° - 50°
- Adjustable 6" shear reduction recline
- Powered by MK batteries, rated up to 19 miles per charge

Accessories

(Based on chair configuration)

- USB charger
- Oxygen tank holder
- Cup holder

Options

- ASL alternative drive
- Therafin positioning products
- Power articulating footplate
- Push handles
- Proform headrest with EZ-Clamp bracket



Features MK batteries



Specifications

Product#	P335
Weight Capacity	300 lbs.
Drive train	FWD
Seat Width	14-22"
Seat Depth	14-22"
Backrest Height	18-22"
Seat-to-Floor Height	17.5-19.5"
Seat Elevation (Optional)	12"
Overall Length	39.6"
Overall Width	24.5"
Overall Height	41.5"
Turning Radius	26" w/o footplate
Ground Clearance	3"
Total Weight	350-385 lbs.
Max Speed	6 mph
Battery Size (2 required)	MK 60Ah, M34 Gel
Range Per Charge	up to 19 mi.
Motor	2-Pole/222W
Controller	R-NET 120A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	14" x 3" hybrid gel
Rear Wheels	9" hybrid gel
HCPCS Code	K0856/K0861

Color Options

