



AXCEL NARROW ORDER FORM

| Standard Features | Base | Power Tilt | Power Tilt & AFP | Power Tilt, Lift & AFP | Power Tilt, Recline & Lift |
|---|---------------|----------------|------------------|------------------------|----------------------------|
| 14" Drive Wheels | ✓ | ✓ | ✓ | ✓ | ✓ |
| 6" Castors | ✓ | ✓ | ✓ | ✓ | ✓ |
| Choice of Coloured Shroud (Grey, Black, Red, Matte Blue, Green) | ✓ | ✓ | ✓ | ✓ | ✓ |
| Independent, Hydraulic Suspension | ✓ | ✓ | ✓ | ✓ | ✓ |
| 222W Motors | ✓ | ✓ | ✓ | ✓ | ✓ |
| MK34 60Ah Batteries | ✓ | ✓ | ✓ | ✓ | ✓ |
| Centre Footplate & Mount | ✓ | ✓ | ✓ | ✓ | ✓ |
| Swing Away Joystick Bracket | ✓ | ✓ | ✓ | ✓ | ✓ |
| JSM (LED-L) Joystick | ✓ | ✓ | ✓ | ✓ | ✓ |
| Gyro Stabilized Ride Control | ✓ | ✓ | ✓ | ✓ | ✓ |
| Multi Axis Headrest, Pad and Bracket | Opt | ✓ | ✓ | ✓ | ✓ |
| Manual Tilt & Recline Seat | ✓ | N/A | N/A | N/A | N/A |
| Frequency Power Tilt only Seat | N/A | ✓ | N/A | ✓ | N/A |
| Frequency Power Tilt & Recline Anti-Shear Seat | N/A | N/A | ✓ | N/A | ✓ |
| 12" Powered Seat Lift | N/A | N/A | N/A | ✓ | ✓ |
| T-Bar Adjustable Armrests | ✓ | N/A | N/A | N/A | N/A |
| Cantilever Folding Armrests | N/A | ✓ | ✓ | ✓ | ✓ |
| Seat Belt | ✓ | ✓ | ✓ | ✓ | ✓ |
| LED Front and Rear Lights | ✓ | ✓ | ✓ | ✓ | ✓ |
| Merits Seat Cushion | ✓ | N/A | N/A | N/A | N/A |
| Comfort Company Curve Seat Cushion | N/A | ✓ | ✓ | ✓ | ✓ |
| Merits Seat Back | ✓ | N/A | N/A | N/A | N/A |
| Comfort Company Acta Back | N/A | ✓ | ✓ | ✓ | ✓ |
| Price | £9,000 | £11,000 | £12,200 | £12,000 | £13,200 |

AXCEL NARROW ORDER FORM



Account #: _____ Date: _____

Dealer Name: _____

Contact: _____

Phone _____

Email: _____


PO Number: _____

Ship to Address: _____

City: _____

Region /County: _____

Post Code: _____



| | | Seat/Leg Length | | | | | | | | | |
|------|--------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--|
| | | Wide | | | | | | | | | |
| | Inches | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | |
| Deep | 14 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |
| | 15 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |
| | 16 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |
| | 17 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |
| | 18 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |
| | 19 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |
| | 20 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |
| | 21 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |
| | 22 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | |

CLIENT INFORMATION





- A. Hip Width _____
- B. Seat Depth _____
- C. Lower Leg Length _____
- D. Seat to Top of Head _____
- E. Seat to Top of Shoulder _____
- F. Seat to Axilla (L, R) Armpit _____
- G. Seat to Arm _____
- H. Chest Width _____
- I. Shoulder Width _____
- J. Client Height _____
- K. Client Weight _____

| Availability | Current Name | Width | Depth | Phys Length |
|--------------------------------|--------------|----------|-----------------|-------------|
| <input type="checkbox"/> Dev F | N/A | 14 - 15" | 14" Deep | 17" long |
| <input type="checkbox"/> Dev F | N/A | 14 - 15" | 15 - 14.5" Deep | 18" long |
| <input type="checkbox"/> Avail | Short | 14 - 15" | 16 - 15.5" Deep | 18.25" Long |
| <input type="checkbox"/> Avail | Long | 16 - 22" | 16 - 22" Deep | 24" long |
| Configuration Not Available | | | | N/A |

| Seating Adjustments - REQUIRED | | | |
|--------------------------------------|------------------------------|------------------------------|-------------|
| Seat Back Width | Seat Pan Width | Depth Pan Seat Rail | |
| <input type="checkbox"/> 16" | <input type="checkbox"/> 14" | <input type="checkbox"/> 14" | 42610013 |
| <input type="checkbox"/> 18" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | 42610016-NI |
| <input type="checkbox"/> 20" | <input type="checkbox"/> 16" | <input type="checkbox"/> 18" | 42610016-NI |
| <input type="checkbox"/> 22" | <input type="checkbox"/> 17" | <input type="checkbox"/> 20" | 42610016-NI |
| <input type="checkbox"/> Aftermarket | <input type="checkbox"/> 18" | <input type="checkbox"/> 22" | 42610016-NI |
| | <input type="checkbox"/> 19" | | |
| | <input type="checkbox"/> 20" | | |
| | <input type="checkbox"/> 21" | | |
| | <input type="checkbox"/> 22" | | |

| Angle-Adjustable Back - If not selected, it will be set to 95° | | | | |
|--|--------------------------------|--------------------------------|-------------------------------|-------------------------------|
| <input type="checkbox"/> 90° | <input type="checkbox"/> 95° | <input type="checkbox"/> 100° | <input type="checkbox"/> 105° | <input type="checkbox"/> 110° |
| Only applicable for power tilt only. | | | | |
| Seat to Floor Height Options - REQUIRED | | | | |
| <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18.5" | <input type="checkbox"/> 19.5" | | |

CUSTOM OPTIONS

| Color Options - REQUIRED | | Part # | £ MSRP |
|--------------------------|---|----------|--------|
| <input type="checkbox"/> |  Ruby Red P335 | 30500970 | £0.00 |
| <input type="checkbox"/> |  Noir Black P335 | 30500969 | £0.00 |
| <input type="checkbox"/> |  Camo Green P335 | 31300010 | £0.00 |
| <input type="checkbox"/> |  Dove Grey P335 | 30500951 | £0.00 |
| <input type="checkbox"/> |  Sapphire Blue, P335 | 31300009 | £0.00 |

| Batteries - REQUIRED | | Part # | £ MSRP |
|--------------------------|---------------------------------|--------------|---------|
| <input type="checkbox"/> | Battery, M34, 60Ah, (Gel) Rehab | M50-12 SLD M | £354.33 |

| Joystick Type - REQUIRED | | Part # | £ MSRP |
|--------------------------|--|--------|-----------|
| <input type="checkbox"/> | Joystick Cjsm2-LW/With Lighting, R-Net | D51595 | £940.94 |
| <input type="checkbox"/> | CJSM2 w/Bluetooth | D51634 | £1,019.69 |

| Expansion Modules | | Part # | £ MSRP |
|--------------------------|---|------------|---------|
| <input type="checkbox"/> | R-Net Omni 2 Display Section & Cable | D51709 | £708.66 |
| <input type="checkbox"/> | R-Net Omni 2 Input Section | D51710 | £629.92 |
| <input type="checkbox"/> | R-Net Omni 2 Cable 1.5M | SA81774L15 | £54.33 |
| <input type="checkbox"/> | Swing Away Joystick Display Mount, Omni | ARMS260-Q | £295.28 |
| <input type="checkbox"/> | Display Adapter, Omni 2 | P81996 | £35.34 |

| Attendant Kit | | Part # | £ MSRP |
|--------------------------|-----------------------------------|--------|---------|
| <input type="checkbox"/> | Dual Attendant Module, R-Net P323 | D50882 | £692.13 |
| <input type="checkbox"/> | Bracket, Attendant Control | J3D | £3.46 |

| Optional Programmers | | Part # | £ MSRP |
|--------------------------|----------------------|--------|---------|
| <input type="checkbox"/> | PC Programmer, R-Net | D50612 | £687.31 |

| Joystick Mounting Brackets - REQUIRED | | Part # | £ MSRP |
|---------------------------------------|---|-------------|---------|
| <input type="checkbox"/> | Swing-away Joystick Bracket, Rh (Lt Rehand) | 39922019-NI | £216.54 |
| <input type="checkbox"/> | Swing-away Joystick Bracket, Lh, P323, P325 - Can B | 39922020-NI | £216.54 |
| <input type="checkbox"/> | Bracket Kit, In Line P325 | 33500883 | £17.32 |

| Joystick Knob - OPTIONAL | | Part # | £ MSRP |
|--------------------------|---|---------|---------|
| <input type="checkbox"/> | Joystick Knob, Mushroom (M48) P&G/Vr2 Only - P325 | SA72439 | £177.51 |
| <input type="checkbox"/> | Joystick Knob, 38Mm Ball (B38) | SA72441 | £119.69 |
| <input type="checkbox"/> | Joystick Knob, Soft Ball (Sb70) | SA81095 | £86.61 |

| Armrests | | Part # | £ MSRP |
|--------------------------|---|--------|---------|
| <input type="checkbox"/> | Armrests 2X12 (pair) | 212G | £106.30 |
| <input type="checkbox"/> | Armrests 2X14 (pair) | 214G | £112.60 |
| <input type="checkbox"/> | Armrests 3.5X12 (pair) | 312G | £121.26 |
| <input type="checkbox"/> | Armrests 3.5X14 (pair) | 314G | £147.24 |
| <input type="checkbox"/> | Bracket, 3" X 4", Gel Pad | UESF1 | £92.91 |
| <input type="checkbox"/> | Arm Support, Polyurethane Trough 5" W x 11" L | PUT11 | £133.86 |
| <input type="checkbox"/> | Hand rest Pad, Flat 4.5" W x 7" L for PUT11 | PUTFH | £102.36 |
| <input type="checkbox"/> | Handle Right for PUT11 8.5"L, Anatomic | PUTAHR | £118.11 |
| <input type="checkbox"/> | Handle Left for PUT11 8.5"L, Anatomic | PUTAHL | £118.11 |
| <input type="checkbox"/> | Hand Rest with Palm Extensor for PUT11 | PUHPPE | £118.11 |

| Headrests - REQUIRED | | Part # | £ MSRP |
|--------------------------|---|----------|---------|
| <input type="checkbox"/> | Headrest, Midline Hardware, ProForm | 32100388 | £157.48 |
| <input type="checkbox"/> | Headrest, Multi-Axis Hardware, ProForm (No Pad) | 32100372 | £173.23 |
| <input type="checkbox"/> | Headrest Pad 10, PC Prom | 32000640 | £165.35 |
| <input type="checkbox"/> | Headrest Pad, 14" | 50367 | £123.62 |

| Footplates | | Part # | £ MSRP |
|--------------------------|---|----------|---------|
| <input type="checkbox"/> | Footplate Assembly W/Bracket, Alum (W/ Merits Logo) | 33300350 | £0.00 |
| <input type="checkbox"/> | Center Footplate and Mount-P335 V3 | 33300502 | £75.59 |
| <input type="checkbox"/> | Split Footplate, 5"x6", Pair | 33300481 | £125.98 |
| <input type="checkbox"/> | Split Footplates 3 Inch Extension Kit, Pair | 33300489 | £43.31 |
| <input type="checkbox"/> | Side Stop Plate for 5"x6" Split Footplate, Pair | 33300495 | £22.05 |
| <input type="checkbox"/> | Split Footplate, 5"x11", Pair | 33300473 | £138.58 |
| <input type="checkbox"/> | Split Footplates 3 Inch Extension Kit, Pair | 33300489 | £43.31 |
| <input type="checkbox"/> | Side Stop Plate for 5"x11" Split Footplate, Pair | 33300493 | £22.05 |
| <input type="checkbox"/> | Split Footplate Linkage Bar | 33300568 | £4.72 |

| Lower Extremity Support - OPTIONAL | Part # | £ MSRP |
|--|-----------|---------|
| <input type="checkbox"/> Adductor Bracket Swing-Away Mounting Hardware | SRASF1-NI | £195.28 |
| <input type="checkbox"/> Adductor Knee or Thigh Gel Pad, 3x4, Sta-Rite | SRASF1 | £238.58 |
| <input type="checkbox"/> Abductor Pad, M, P323 | 39900174 | £64.57 |
| <input type="checkbox"/> Abductor Pad, S, P323 | 39900175 | £62.99 |
| <input type="checkbox"/> Swing Away Abductor Bracket, P323 & P335 | 39900176 | £129.13 |
| <input type="checkbox"/> Residual Limb Support, 6x7.5, Sta-Rite | SRSSF | £303.15 |
| <input type="checkbox"/> Residual Limb Support Mounting Hardware, Swing-Away | SRSSFHO | £259.84 |
| <input type="checkbox"/> Knee Adductor, Flip-down WShort Arm, Adapter PR | 30240 | £480.31 |
| <input type="checkbox"/> Knee Adductor, Flip-down WLong Arm, Adapter, PR | 30276 | £480.31 |
| <input type="checkbox"/> Pad, Med Knee Adductor W Neoprene Cover | 50491 | £93.70 |
| <input type="checkbox"/> Pad, SM Knee Adductor W Neoprene Cover | 50490 | £93.70 |
| <input type="checkbox"/> Footsure Ankle Support W Buckle, Med PR | 60252 | £141.73 |
| <input type="checkbox"/> Footsure Ankle Support W Buckle, LG PR | 60253 | £141.73 |
| <input type="checkbox"/> Footsure Ankle Support W Hook-N-Loop, Med Pr. | 60257 | £141.73 |
| <input type="checkbox"/> Footsure Ankle Support W Hook-N-Loop, LG PR | 60258 | £141.73 |
| <input type="checkbox"/> Footbox, Med | 30112 | £252.76 |
| <input type="checkbox"/> Footbox, LG | 30113 | £252.76 |
| <input type="checkbox"/> Footbox, Split Med | 30153 | £318.11 |
| <input type="checkbox"/> Footbox, Split LG | 30154 | £323.62 |
| <input type="checkbox"/> Shoe holder W Plain Straps, Med PR | 30303 | £122.05 |
| <input type="checkbox"/> Shoe holder W Plain Straps, LG PR | 30309 | £122.05 |
| <input type="checkbox"/> Shoe holder W Plain Straps, X-LG PR | 30316 | £122.05 |
| <input type="checkbox"/> Flexsure Feet Padded Inserts, Med, PR | 31762 | £70.87 |
| <input type="checkbox"/> Flexsure Feet Padded Inserts, LG, PR | 31763 | £70.87 |
| <input type="checkbox"/> Flexsure Feet Padded Inserts, X-LG, PR | 31764 | £70.87 |
| <input type="checkbox"/> Flexsure Feet W Hook-N-Loop Closure, Med, PR | 31766 | £252.76 |
| <input type="checkbox"/> Flexsure Feet W Hook-N-Loop Closure, LG, PR | 31767 | £252.76 |
| <input type="checkbox"/> Flexsure Feet W Hook-N-Loop Closure, X-LG, PR | 31768 | £252.76 |
| <input type="checkbox"/> Flexsure Feet W Buckle Closure, Med, PR | 31772 | £252.76 |
| <input type="checkbox"/> Flexsure Feet W Buckle Closure, LG, PR | 31773 | £252.76 |
| <input type="checkbox"/> Flexsure Feet W Buckle Closure, X-LG, PR | 31774 | £252.76 |

| Pelvic Belts - REQUIRED | Part # | £ MSRP |
|---|----------|---------|
| <input type="checkbox"/> Seat Belt 66 Inch (Loop Type) | 43700319 | £31.50 |
| <input type="checkbox"/> Seat Belt W/ Axle Mount 84" | SB225 | £69.29 |
| <input type="checkbox"/> Belt, Therafit 4Pt Y-Style W Buckle, SM | 60024 | £133.07 |
| <input type="checkbox"/> Belt, Therafit 4Pt Y-Style W Buckle, Med | 60025 | £133.07 |
| <input type="checkbox"/> Belt, Therafit 4Pt Y-Style W Buckle, LG | 60026 | £133.07 |
| <input type="checkbox"/> Belt, 1.5 W Single Pull, 5.25W X 2.25H Pads | 62001 | £71.65 |
| <input type="checkbox"/> Belt, 1.5 W Single Pull, 7.25W X 2.25H Pads | 62101 | £71.65 |
| <input type="checkbox"/> Belt, 2 W Single Pull, 9.25W X 3H Pads W Clips | 64001 | £76.38 |

| Pelvic Belts - REQUIRED | | Part # | £ MSRP |
|--------------------------|---|--------|--------|
| <input type="checkbox"/> | Belt, 2 W Single Pull, 11.25W X 3H Pads W Clips | 64101 | £76.38 |

| Upper Extremity Support - OPTIONAL | | Part # | £ MSRP |
|------------------------------------|--|--------|---------|
| <input type="checkbox"/> | Lateral Swing-away kit (including 0 & 1 Offset), T | 50560 | £211.02 |
| <input type="checkbox"/> | Lateral Pad, Flat 3 x 5, Therafin | 50052 | £66.93 |
| <input type="checkbox"/> | Lts Pad, 4 X 5, Flat | 50053 | £66.93 |
| <input type="checkbox"/> | Lateral Pad, Flat 4x6, Therafin | 50054 | £66.93 |
| <input type="checkbox"/> | Lateral Pad, Curved 3x5, Therafin | 50057 | £83.46 |
| <input type="checkbox"/> | Lts Pad, 4 X 5, Curved | 50058 | £83.46 |
| <input type="checkbox"/> | Lateral Pad, Curved 4x6, Therafin | 50059 | £83.46 |
| <input type="checkbox"/> | Vest, Comfort fit, Full/Med | 30590 | £126.77 |
| <input type="checkbox"/> | Vest, Comfort fit, Full/LG | 30591 | £126.77 |
| <input type="checkbox"/> | Vest, Comfort fit, Full/X-LG | 30592 | £126.77 |
| <input type="checkbox"/> | Vest, Comfort fit, Trim/Med | 30619 | £128.35 |
| <input type="checkbox"/> | Vest, Comfort fit, Trim/LG | 30620 | £128.35 |
| <input type="checkbox"/> | Vest, Comfort fit, Trim/X-LG | 30621 | £133.07 |
| <input type="checkbox"/> | Shoulder Harness, Therafit, Med | 61503 | £123.62 |
| <input type="checkbox"/> | Shoulder Harness, Therafit, LG | 61703 | £129.13 |
| <input type="checkbox"/> | Shoulder Harness, Therafit, X-LG | 61803 | £134.65 |
| <input type="checkbox"/> | Chest Strap, Therafit Dynamic W Single D-Ring Med | 60114 | £105.51 |
| <input type="checkbox"/> | Chest Strap, Therafit Dynamic W Single D-Ring LG | 60115 | £105.51 |
| <input type="checkbox"/> | Chest Strap, Therafit Dynamic W Single D-Ring X-LG | 60123 | £122.83 |
| <input type="checkbox"/> | Chest Strap, Therafit Dynamic W Hook-N-Loop Med | 60120 | £105.51 |
| <input type="checkbox"/> | Chest Strap, Therafit Dynamic W Hook-N-Loop LG | 60121 | £105.51 |
| <input type="checkbox"/> | Strap, Positioning, Universal, 4.5W X 56L, Med | 34435 | £55.12 |
| <input type="checkbox"/> | Strap, Positioning, Universal, 6W X 56L LG | 34436 | £59.06 |

| Accessories | | Part # | £ MSRP |
|--------------------------|---|--------------|-----------|
| <input type="checkbox"/> | Strap, Positioning, Universal, 6W X 56L LG | 34436 | £59.06 |
| <input type="checkbox"/> | Micro Extremity Control - Joystick Only | ASL 130 | £2,449.61 |
| <input type="checkbox"/> | Nucleus Mount | ASL 660 | £543.31 |
| <input type="checkbox"/> | Joystick Mount Adj. Height Platform For Mec | ASL 603M | £543.31 |
| <input type="checkbox"/> | Joystick Mount Aluminum Bib (Standard) | ASL 625 | £748.03 |
| <input type="checkbox"/> | Midline Palm Tray With Mec Insert Aluminum | ASL 671 | £748.03 |
| <input type="checkbox"/> | Headrest Mount Swing-Away Adjustable | ASL 628 | £543.31 |
| <input type="checkbox"/> | Headrest Mounting Hardware (Adj/Rem)P323 | ASL 606 | £303.15 |
| <input type="checkbox"/> | Headrest 10" , No Mounting Hardware | ASL 920 | £354.33 |
| <input type="checkbox"/> | Head Array ASL Atom Electronic , Fixed Mount Hdwr | ASL 104 ATOM | £3,304.96 |
| <input type="checkbox"/> | Headrest Mounting Hardware (Adj/Rem)P323 | ASL 606 | £303.15 |
| <input type="checkbox"/> | Headrest Mounting Hardware (Adj/Rem)P323 | ASL 606 | £303.15 |
| <input type="checkbox"/> | Digital Head Array ASL Fusion | ASL 110 | £4,129.76 |

| Accessories | Part # | £ MSRP |
|--|-------------|---------|
| <input type="checkbox"/> Occipital Pad For The ASL 110 Fusion | ASL 935 | £108.66 |
| <input type="checkbox"/> Swing-Away Hardware For Lateral Pad 9 Inch LH | ASL 932UL-S | £174.80 |
| <input type="checkbox"/> Lateral Spot Pad With Vinyl Cover - Lh | ASL 933VL | £108.66 |
| <input type="checkbox"/> Swing-Away Hardware For Lateral Pad 9 Inch RH | ASL 932UR-S | £174.80 |
| <input type="checkbox"/> Lateral Spot Pad With Vinyl Cover - Rh | ASL 933VR | £108.66 |
| <input type="checkbox"/> Headrest Mounting Hardware (Adj/Rem)P323 | ASL 606 | £303.15 |
| <input type="checkbox"/> Fusion Display And Programmer | ASL 165 | £347.24 |
| <input type="checkbox"/> Pulse Mechanical Switch - Black | ASL 320 BK | £122.83 |
| <input type="checkbox"/> Electronic Proximity Sensor Switch - 3-Pin | ASL 204 | £519.69 |
| <input type="checkbox"/> Headrest Swing-away,Adjustable Switch Mount Assy Rh | ASL 930L | £382.68 |
| <input type="checkbox"/> Headrest Swing-away,Adjustable Switch Mount Assy Lh | ASL 930R | £382.68 |
| <input type="checkbox"/> Tray, Clear, 25.5W X 22D W 15X9 Rim & Clamps | T16R12A | £342.52 |
| <input type="checkbox"/> Tray, Clear, 27.5W X 23.5D W 16X10 Rim & Clamps | T17L12A | £362.20 |
| <input type="checkbox"/> Tray, Clear, 27.5W X 23.5D W 16X10Rim & Clamps | T17R12A | £362.20 |
| <input type="checkbox"/> Half-Tray, Padded Trimline , Flip LH | 41351 | £212.60 |
| <input type="checkbox"/> Half-Tray, Padded Trimline , Flip RH | 41352 | £212.60 |
| <input type="checkbox"/> Half-Tray, Padded Trimline W/Cupholder, Flip RH | 41307B | £227.56 |
| <input type="checkbox"/> Half-Tray, Padded Trimline W/Cupholder, Flip LH | 41308B | £227.56 |
| <input type="checkbox"/> Half-Tray, Viewline+, Flip LH | 41311 | £193.70 |
| <input type="checkbox"/> Half-Tray, Viewline+, Flip RH | 41310 | £193.70 |
| <input type="checkbox"/> Oxygen Holder P335 | 34100319 | £173.23 |
| <input type="checkbox"/> Cup Holder-Manual WC - Pow/Scoot | 46700060 | £20.47 |
| <input type="checkbox"/> Universal Cup Holder Bracket, Gen 2 | 34100323 | £31.49 |
| <input type="checkbox"/> USB Charger - Universal | 61230167 | £35.43 |
| <input type="checkbox"/> Push Handle -P335 (Set) | 34100328 | £173.23 |

| Comfort Company Curve Seat Cushion | Part # | £ MSRP |
|--|---------------|---------|
| <input type="checkbox"/> Curve Cushion, 14" X 14", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-14W14D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 14" X 15", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-14W15D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 14" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-14W16D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 14" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-14W17D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 14" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-14W18D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 15" X 14", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-15W14D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 15" X 15", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-15W15D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 15" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-15W16D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 15" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-15W17D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 15" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-15W18D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 15" X 19", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-15W19D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 15" X 20", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-15W20D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 16" X 14", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W14D | £100.00 |
| <input type="checkbox"/> Curve Cushion, 16" X 15", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W15D | £100.00 |

| Comfort Company Curve Seat Cushion | | Part # | £ MSRP |
|------------------------------------|---|---------------|---------|
| <input type="checkbox"/> | Curve Cushion, 16" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W16D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 16" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W17D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 16" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W18D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 16" X 19", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W19D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 16" X 20", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W20D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 16" X 21", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W21D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 16" X 22", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-16W22D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 17" X 15", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-17W15D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 17" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-17W16D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 17" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-17W17D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 17" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-17W18D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 17" X 19", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-17W19D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 17" X 20", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-17W20D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 18" X 15", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-18W15D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 18" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-18W16D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 18" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-18W17D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 18" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-18W18D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 18" X 19", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-18W19D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 18" X 20", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-18W20D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 18" X 21", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-18W21D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 18" X 22", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-18W22D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 19" X 15", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-19W16D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 19" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-19W17D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 19" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-19W18D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 19" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-19W19D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 19" X 19", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-19W20D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 19" X 20", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-20W16D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 19" X 21", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-20W17D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 19" X 22", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-20W18D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 20" X 15", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-20W19D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 20" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-20W20D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 20" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-20W21D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 20" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-20W22D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 20" X 19", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-21W16D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 20" X 20", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-21W17D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 20" X 21", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-21W18D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 21" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-21W19D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 21" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-21W20D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 21" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-22W16D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 21" X 19", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-22W17D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 21" X 20", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-22W18D | £100.00 |

| Comfort Company Curve Seat Cushion | | Part # | £ MSRP |
|------------------------------------|---|---------------|---------|
| <input type="checkbox"/> | Curve Cushion, 21" X 21", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-22W19D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 22" X 16", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-22W20D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 22" X 17", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-22W21D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 22" X 18", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-22W22D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 22" X 19", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-23W16D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 22" X 20", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-23W17D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 22" X 21", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-23W18D | £100.00 |
| <input type="checkbox"/> | Curve Cushion, 22" X 22", Molded Base, FR Stretch-Air Cover w/ Non-Skid & Loop Bottom | CU-FRS-23W19D | £100.00 |

| Comfort Company Acta Embrace Seat Cushion | | Part # | £ MSRP |
|---|---|----------------|---------|
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 14" X 14", Fr Stretch-Air | EM-FRSA-14W14D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 14" X 15", Fr Stretch-Air | EM-FRSA-14W15D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 14" X 16", Fr Stretch-Air | EM-FRSA-14W16D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 14" X 17", Fr Stretch-Air | EM-FRSA-14W17D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 14" X 18", Fr Stretch-Air | EM-FRSA-14W18D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 14" X 19", Fr Stretch-Air | EM-FRSA-14W19D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 14" X 20", Fr Stretch-Air | EM-FRSA-14W20D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 10", Fr Stretch-Air | EM-FRSA-15W10D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 11", Fr Stretch-Air | EM-FRSA-15W11D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 12", Fr Stretch-Air | EM-FRSA-15W12D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 13", Fr Stretch-Air | EM-FRSA-15W13D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 14", Fr Stretch-Air | EM-FRSA-15W14D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 15", Fr Stretch-Air | EM-FRSA-15W15D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 16", Fr Stretch-Air | EM-FRSA-15W16D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 17", Fr Stretch-Air | EM-FRSA-15W17D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 18", Fr Stretch-Air | EM-FRSA-15W18D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 19", Fr Stretch-Air | EM-FRSA-15W19D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 20", Fr Stretch-Air | EM-FRSA-15W20D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 15" X 21", Fr Stretch-Air | EM-FRSA-15W21D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 12", Fr Stretch-Air | EM-FRSA-16W12D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 13", Fr Stretch-Air | EM-FRSA-16W13D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 14", Fr Stretch-Air | EM-FRSA-16W14D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 15", Fr Stretch-Air | EM-FRSA-16W15D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 16", Fr Stretch-Air | EM-FRSA-16W16D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 17", Fr Stretch-Air | EM-FRSA-16W17D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 18", Fr Stretch-Air | EM-FRSA-16W18D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 19", Fr Stretch-Air | EM-FRSA-16W19D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 20", Fr Stretch-Air | EM-FRSA-16W20D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 21", Fr Stretch-Air | EM-FRSA-16W21D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 16" X 22", Fr Stretch-Air | EM-FRSA-16W22D | £313.00 |
| <input type="checkbox"/> | Cover: Embrace Anti-Thrust Cushion, 17" X 12", Fr Stretch-Air | EM-FRSA-17W12D | £313.00 |

| Comfort Company Acta Embrace Seat Cushion | Part # | £ MSRP |
|--|----------------|---------|
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 13", Fr Stretch-Air | EM-FRSA-17W13D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 14", Fr Stretch-Air | EM-FRSA-17W14D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 15", Fr Stretch-Air | EM-FRSA-17W15D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 16", Fr Stretch-Air | EM-FRSA-17W16D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 17", Fr Stretch-Air | EM-FRSA-17W17D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 18", Fr Stretch-Air | EM-FRSA-17W18D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 19", Fr Stretch-Air | EM-FRSA-17W19D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 20", Fr Stretch-Air | EM-FRSA-17W20D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 21", Fr Stretch-Air | EM-FRSA-17W21D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 17" X 22", Fr Stretch-Air | EM-FRSA-17W22D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 14", Fr Stretch-Air | EM-FRSA-18W14D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 15", Fr Stretch-Air | EM-FRSA-18W15D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 16", Fr Stretch-Air | EM-FRSA-18W16D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 17", Fr Stretch-Air | EM-FRSA-18W17D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 18", Fr Stretch-Air | EM-FRSA-18W18D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 19", Fr Stretch-Air | EM-FRSA-18W19D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 20", Fr Stretch-Air | EM-FRSA-18W20D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 21", Fr Stretch-Air | EM-FRSA-18W21D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 22", Fr Stretch-Air | EM-FRSA-18W22D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 23", Fr Stretch-Air | EM-FRSA-18W23D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 18" X 24", Fr Stretch-Air | EM-FRSA-18W24D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 14", Fr Stretch-Air | EM-FRSA-19W14D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 15", Fr Stretch-Air | EM-FRSA-19W15D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 16", Fr Stretch-Air | EM-FRSA-19W16D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 17", Fr Stretch-Air | EM-FRSA-19W17D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 18", Fr Stretch-Air | EM-FRSA-19W18D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 19", Fr Stretch-Air | EM-FRSA-19W19D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 20", Fr Stretch-Air | EM-FRSA-19W20D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 21", Fr Stretch-Air | EM-FRSA-19W21D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 22", Fr Stretch-Air | EM-FRSA-19W22D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 23", Fr Stretch-Air | EM-FRSA-19W23D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 19" X 24", Fr Stretch-Air | EM-FRSA-19W24D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 15", Fr Stretch-Air | EM-FRSA-20W15D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 16", Fr Stretch-Air | EM-FRSA-20W16D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 17", Fr Stretch-Air | EM-FRSA-20W17D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 18", Fr Stretch-Air | EM-FRSA-20W18D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 19", Fr Stretch-Air | EM-FRSA-20W19D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 20", Fr Stretch-Air | EM-FRSA-20W20D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 21", Fr Stretch-Air | EM-FRSA-20W21D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 22", Fr Stretch-Air | EM-FRSA-20W22D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 23", Fr Stretch-Air | EM-FRSA-20W23D | £313.00 |

| Comfort Company Acta Embrace Seat Cushion | Part # | £ MSRP |
|--|----------------|---------|
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 20" X 24", Fr Stretch-Air | EM-FRFA-20W24D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 16", Fr Stretch-Air | EM-FRSA-21W16D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 17", Fr Stretch-Air | EM-FRSA-21W17D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 18", Fr Stretch-Air | EM-FRSA-21W18D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 19", Fr Stretch-Air | EM-FRSA-21W19D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 20", Fr Stretch-Air | EM-FRSA-21W20D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 21", Fr Stretch-Air | EM-FRSA-21W21D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 22", Fr Stretch-Air | EM-FRSA-21W22D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 23", Fr Stretch-Air | EM-FRSA-21W23D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 21" X 24", Fr Stretch-Air | EM-FRSA-21W24D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 22" X 16", Fr Stretch-Air | EM-FRSA-22W16D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 22" X 17", Fr Stretch-Air | EM-FRSA-22W17D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 22" X 18", Fr Stretch-Air | EM-FRSA-22W18D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 22" X 19", Fr Stretch-Air | EM-FRSA-22W19D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 22" X 20", Fr Stretch-Air | EM-FRSA-22W20D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 22" X 21", Fr Stretch-Air | EM-FRFA-22W21D | £313.00 |
| <input type="checkbox"/> Cover: Embrace Anti-Thrust Cushion, 22" X 22", Fr Stretch-Air | EM-FRSA-22W22D | £313.00 |

| Comfort Company M2 Seat Cushion | Part # | £ MSRP |
|---|----------------|--------|
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 14" X 14", Fr Stretch-Air | M2-FRSA-14W14D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 14" X 15", Fr Stretch-Air | M2-FRSA-14W15D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 14" X 16", Fr Stretch-Air | M2-FRSA-14W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 14" X 17", Fr Stretch-Air | M2-FRSA-14W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 14" X 18", Fr Stretch-Air | M2-FRSA-14W18D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 14" X 19", Fr Stretch-Air | M2-FRSA-14W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 14" X 20", Fr Stretch-Air | M2-FRSA-14W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 10", Fr Stretch-Air | M2-FRSA-15W10D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 11", Fr Stretch-Air | M2-FRSA-15W11D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 12", Fr Stretch-Air | M2-FRSA-15W12D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 13", Fr Stretch-Air | M2-FRSA-15W13D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 14", Fr Stretch-Air | M2-FRSA-15W14D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 15", Fr Stretch-Air | M2-FRSA-15W15D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 16", Fr Stretch-Air | M2-FRSA-15W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 17", Fr Stretch-Air | M2-FRSA-15W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 18", Fr Stretch-Air | M2-FRSA-15W18D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 19", Fr Stretch-Air | M2-FRSA-15W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 20", Fr Stretch-Air | M2-FRSA-15W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 15" X 21", Fr Stretch-Air | M2-FRSA-15W21D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 12", Fr Stretch-Air | M2-FRSA-16W12D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 13", Fr Stretch-Air | M2-FRSA-16W13D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 14", Fr Stretch-Air | M2-FRSA-16W14D | £392 |

| Comfort Company M2 Seat Cushion | Part # | £ MSRP |
|---|-------------------|--------|
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 15", Fr Stretch-Air | M2-FRSA-16W15D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 16", Fr Stretch-Air | M2-FRSA-16W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 17", Fr Stretch-Air | M2-FRSA-16W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 18", Fr Stretch-Air | M2-FRSA-16W18D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 19", Fr Stretch-Air | M2-FRSA-16W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 20", Fr Stretch-Air | M2-FRSA-16W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 21", Fr Stretch-Air | M2-FRSA-16W21D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 16" X 22", Fr Stretch-Air | M2-FRSA-16W22D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 12", Fr Stretch-Air | M2-FRSA-17W12D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 13", Fr Stretch-Air | M2-FRSA-17W13D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 14", Fr Stretch-Air | M2-FRSA-17W14D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 15", Fr Stretch-Air | M2-FRSA-17W15D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 16", Fr Stretch-Air | M2-FRSA-17W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 17", Fr Stretch-Air | M2-FRSA-17W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 18", Fr Stretch-Air | M2-FRSA-17W18D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 19", Fr Stretch-Air | M2-FRSA-17W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 20", Fr Stretch-Air | M2-FRSA-17W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 21", Fr Stretch-Air | M2-FRSA-17W21D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 17" X 22", Fr Stretch-Air | M2-FRSA-17W22D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 14", Fr Stretch-Air | M2-FRSA-18W14D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 15", Fr Stretch-Air | M2-FRSA-18W15D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 16", Fr Stretch-Air | M2-FRSA-18W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 17", Fr Stretch-Air | M2-FRSA-18W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 18", Fr Stretch-Air | M2-FRSA-18W18D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 19", Fr Stretch-Air | M2-FRSA-18W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 20", Fr Stretch-Air | M2-FRSA-18W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 21", Fr Stretch-Air | M2-FRSA-18W21D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 22", Fr Stretch-Air | M2-FRSA-18W22D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 23", Fr Stretch-Air | M2-FRSA-18W23D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 18" X 24", Fr Stretch-Air | M2-FRSA-18W24D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 14", Fr Stretch-Air | M2-FRSA-19W14D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 15", Fr Stretch-Air | M2-FRSA-19W15D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 16", Fr Stretch-Air | M2-FRSA-19W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 17", Fr Stretch-Air | M2-FRSA-19W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 18", Fr Stretch-Air | M2-FRSA-19W18D-CV | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 19", Fr Stretch-Air | M2-FRSA-19W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 20", Fr Stretch-Air | M2-FRSA-19W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 21", Fr Stretch-Air | M2-FRSA-19W21D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 22", Fr Stretch-Air | M2-FRSA-19W22D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 23", Fr Stretch-Air | M2-FRSA-19W23D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 19" X 24", Fr Stretch-Air | M2-FRSA-19W24D | £392 |

| Comfort Company M2 Seat Cushion | Part # | £ MSRP |
|---|----------------|--------|
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 15", Fr Stretch-Air | M2-FRSA-20W15D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 16", Fr Stretch-Air | M2-FRSA-20W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 17", Fr Stretch-Air | M2-FRSA-20W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 18", Fr Stretch-Air | M2-FRSA-20W18D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 19", Fr Stretch-Air | M2-FRSA-20W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 20", Fr Stretch-Air | M2-FRSA-20W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 21", Fr Stretch-Air | M2-FRSA-20W21D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 22", Fr Stretch-Air | M2-FRSA-20W22D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 23", Fr Stretch-Air | M2-FRSA-20W23D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 20" X 24", Fr Stretch-Air | M2-FRSA-20W24D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 16", Fr Stretch-Air | M2-FRSA-21W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 17", Fr Stretch-Air | M2-FRSA-21W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 18", Fr Stretch-Air | M2-FRSA-21W18D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 19", Fr Stretch-Air | M2-FRSA-21W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 20", Fr Stretch-Air | M2-FRSA-21W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 21", Fr Stretch-Air | M2-FRSA-21W21D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 22", Fr Stretch-Air | M2-FRSA-21W22D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 23", Fr Stretch-Air | M2-FRSA-21W23D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 21" X 24", Fr Stretch-Air | M2-FRSA-21W24D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 22" X 16", Fr Stretch-Air | M2-FRSA-22W16D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 22" X 17", Fr Stretch-Air | M2-FRSA-22W17D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 22" X 18", Fr Stretch-Air | M2-FRSA-22W18D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 22" X 19", Fr Stretch-Air | M2-FRSA-22W19D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 22" X 20", Fr Stretch-Air | M2-FRSA-22W20D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 22" X 21", Fr Stretch-Air | M2-FRSA-22W21D | £392 |
| <input type="checkbox"/> Cover: M2 Anti-Thrust Cushion, 22" X 22", Fr Stretch-Air | M2-FRSA-22W22D | £392 |

| Comfort Company Acta Embrace Seat Back and Cushion | Part # | £ MSRP |
|---|--------------|--------|
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 14W X 16"L, For EMEA | AB-SP-14W16L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 14W X 18"L, For EMEA | AB-SP-14W18L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 16W X 16"L, For EMEA | AB-SP-16W16L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 16W X 18"L, For EMEA | AB-SP-16W18L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 16W X 20"L, For EMEA | AB-SP-16W20L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 16W X 22"L, For EMEA | AB-SP-16W22L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 18W X 16"L, For EMEA | AB-SP-18W16L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 18W X 18"L, For EMEA | AB-SP-18W18L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 18W X 20"L, For EMEA | AB-SP-18W20L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 18W X 22"L, For EMEA | AB-SP-18W22L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 20W X 18"L, For EMEA | AB-SP-20W18L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 20W X 20"L, For EMEA | AB-SP-20W20L | £875 |
| <input type="checkbox"/> Acta-Back: Stretch-Air, Compass Power Mount Corpus, 20W X 22"L, For EMEA | AB-SP-20W22L | £875 |

| Comfort Company Acta Embrace Seat Back and Cushion | | Part # | £ MSRP |
|--|--|--------------|--------|
| <input type="checkbox"/> | Acta-Back: Stretch-Air, Compass Power Mount Corpus, 22W X 18"L, For EMEA | AB-SP-22W18L | £875 |
| <input type="checkbox"/> | Acta-Back: Stretch-Air, Compass Power Mount Corpus, 22W X 20"L, For EMEA | AB-SP-22W20L | £875 |