

Axcel Narrow Specifications

Imperial Measurements

Product #	P335
Weight Capacity	136.08kg
Drivetrain	FWD
Seat Width	355.60mm - 508.00mm
Seat Depth	355.60mm - 508.00mm
Backrest Height	457.20mm - 558.80mm
Seat-to-Floor Height	444.50mm - 495.30mm
Seal Elevation (Opt.)	304.80mm
Overall Length	1005.8mm
Overall Width	571.50mm
Overall Height	1054.1mm
Turning Radius	660.40mm
Ground Clearance	76.200mm
Total Weight	158.76kg - 174.63kg
Max Speed	8.368589kph
Battery Size (2 req.)	MK 60Ah, M34 Gel
Range Per Charge	up to 30.578km
Motor	2 - Pole / 222W
Controller	R-NET 90A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	355.60mm x 76.200mm hybrid gel
Rear Wheels	228.60mm hybrid gel
HCPCS	K0856 / K0861
Standard Colors	Noir Black, Dove Grey, Sapphire Blue, Ruby Red & Camo Green

Metric Measurements

Product #	P335
Weight Capacity	300lb
Drivetrain	FWD
Seat Width	14in - 20in
Seat Depth	14in - 20in
Backrest Height	18in - 22in
Seat-to-Floor Height	17.5in - 19.5in
Seal Elevation (Opt.)	12in
Overall Length	39.6in
Overall Width	22.5in
Overall Height	41.5in
Turning Radius	26in without footplate
Ground Clearance	3in
Total Weight	350lb - 385lb
Max Speed	5.2mph

Battery Size (2 req.)	MK 60Ah, M34 Gel
Range Per Charge	up to 19mi
Motor	2 - Pole / 222W
Controller	R-NET 90A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	14in x 3in hybrid gel
Rear Wheels	9in hybrid gel
HCPCS	K0856 / K0861
Standard Colors	Noir Black, Dove Grey, Sapphire Blue, Ruby Red & Camo Green

