

## Axcel Specifications

### Imperial Measurements

---

<b>Product #</b>	<b>P335</b>
<b>Weight Capacity</b>	136.08kg
<b>Drivetrain</b>	FWD
<b>Seat Width</b>	355.60mm - 558.80mm
<b>Seat Depth</b>	355.60mm - 558.80mm
<b>Backrest Height</b>	457.20mm - 558.80mm
<b>Seat-to-Floor Height</b>	444.50mm - 495.30mm
<b>Seal Elevation (Opt.)</b>	304.80mm
<b>Overall Length</b>	1005.8mm
<b>Overall Width</b>	622.30mm
<b>Overall Height</b>	1054.1mm
<b>Turning Radius</b>	660.40mm
<b>Ground Clearance</b>	76.200mm
<b>Total Weight</b>	158.76kg - 174.63kg
<b>Max Speed</b>	9.656064kph
<b>Battery Size (2 req.)</b>	MK 34 60AH
<b>Range Per Charge</b>	up to 30.578km
<b>Motor</b>	2 - Pole / 222W
<b>Controller</b>	R-NET 120A
<b>Charger</b>	6 Amp
<b>Gradient</b>	7.5°
<b>Front Wheels</b>	Anti-Tipper
<b>Drive Wheels</b>	355.60mm x 76.200mm hybrid gel
<b>Rear Wheels</b>	228.60mm
<b>HCPCS</b>	K0856 / K0861
<b>Standard Colors</b>	Noir Black, Dove Grey, Sapphire Blue, Ruby Red & Camo Green

### Metric Measurements

---

<b>Product #</b>	<b>P335</b>
<b>Weight Capacity</b>	300lb
<b>Drivetrain</b>	FWD
<b>Seat Width</b>	14in - 22in
<b>Seat Depth</b>	14in - 22in
<b>Backrest Height</b>	18in - 22in
<b>Seat-to-Floor Height</b>	17.5in - 19.5in
<b>Seal Elevation (Opt.)</b>	12in
<b>Overall Length</b>	39.6in
<b>Overall Width</b>	24.5in
<b>Overall Height</b>	41.5in
<b>Turning Radius</b>	26in without footplate
<b>Ground Clearance</b>	3in
<b>Total Weight</b>	350lb - 385lb
<b>Max Speed</b>	6mph

<b>Battery Size (2 req.)</b>	MK 34 60AH
<b>Range Per Charge</b>	up to 19mi
<b>Motor</b>	2 - Pole / 222W
<b>Controller</b>	R-NET 120A
<b>Charger</b>	6 Amp
<b>Gradient</b>	7.5°
<b>Front Wheels</b>	Anti-Tipper
<b>Drive Wheels</b>	14in x 3in hybrid gel
<b>Rear Wheels</b>	9in
<b>HCPCS</b>	K0856 / K0861
<b>Standard Colors</b>	Noir Black, Dove Grey, Sapphire Blue, Ruby Red & Camo Green

