

## Ultra Specifications

### Imperial Measurements

---

<b>Product #</b>	<b>P325</b>
<b>Weight Capacity</b>	136.08kg
<b>Drivetrain</b>	MWD
<b>Seat Width</b>	406.40mm - 558.80mm
<b>Seat Depth</b>	406.40mm - 558.80mm
<b>Backrest Height</b>	457.20mm - 711.20mm
<b>Seat-to-Floor Height</b>	Power tilt or pan: 419.10mm - 495.30mm, Captain: 508mm - 584.2mm
<b>Seal Elevation (Opt.)</b>	304.80mm
<b>Tilt (Opt.)</b>	0° - 50°
<b>Recline (Opt.)</b>	90° - 110° Fixed
<b>Articulating Footplate</b>	Optional
<b>Overall Length</b>	914.40mm
<b>Overall Width</b>	609.60mm
<b>Overall Height</b>	1016.0mm
<b>Turning Radius</b>	660.40mm
<b>Ground Clearance</b>	76.200mm
<b>Total Weight</b>	142.88kg - 174.18kg
<b>Max Speed</b>	8.046720kph
<b>Battery Size (2 req.)</b>	34Ah / U1 or 55Ah / 22NF
<b>Range Per Charge</b>	up to 30.578km
<b>Motor</b>	2 - Pole / 240W
<b>Controller</b>	P&G VR2 - A2 60A
<b>Charger</b>	6 Amp
<b>Gradient</b>	6°
<b>Front Wheels</b>	Anti-Tipper
<b>Drive Wheels</b>	355.60mm x 76.200mm
<b>Rear Wheels</b>	152.40mm x 50.800mm
<b>HCPCS</b>	K0822 / K0823 / K0835
<b>Standard Colors</b>	Red, White, Blue & Gray

### Metric Measurements

---

<b>Product #</b>	<b>P325</b>
<b>Weight Capacity</b>	300lb
<b>Drivetrain</b>	MWD
<b>Seat Width</b>	16in - 22in
<b>Seat Depth</b>	16in - 22in
<b>Backrest Height</b>	18in - 28in
<b>Seat-to-Floor Height</b>	Power tilt or pan: 16.5in - 19.5in, Captain: 20in - 23in
<b>Seal Elevation (Opt.)</b>	12in
<b>Tilt (Opt.)</b>	0° - 50°

<b>Recline (Opt.)</b>	90° - 110° Fixed
<b>Articulating Footplate</b>	Optional
<b>Overall Length</b>	36in
<b>Overall Width</b>	24in
<b>Overall Height</b>	40in
<b>Turning Radius</b>	26in
<b>Ground Clearance</b>	2.5in
<b>Total Weight</b>	315lb - 384lb
<b>Max Speed</b>	5mph
<b>Battery Size (2 req.)</b>	34Ah / U1 or 55Ah / 22NF
<b>Range Per Charge</b>	up to 19mi
<b>Motor</b>	2 - Pole / 240W
<b>Controller</b>	P&G VR2 - A2 60A
<b>Charger</b>	6 Amp
<b>Gradient</b>	6°
<b>Front Wheels</b>	Anti-Tipper
<b>Drive Wheels</b>	14in x 3in foam-filled
<b>Rear Wheels</b>	6in x 2in foam-filled
<b>HCPCS</b>	K0822 / K0823 / K0835
<b>Standard Colors</b>	Red, White, Blue & Gray



