

Axcel Narrow Specifications

Imperial Measurements

Product #	P325
Weight Capacity	136.08kg
Drivetrain	MWD
Seat Width	406.40mm - 508.00mm
Seat Depth	406.40mm - 508.00mm
Backrest Height	457.20mm - 711.20mm
Seat-to-Floor Height	419.10mm - 495.30mm
Seal Elevation (Opt.)	304.80mm
Tilt (Opt.)	0° - 50°
Recline (Opt.)	90° - 110° Fixed
Articulating Footplate	Optional
Overall Length	914.40mm
Overall Width	609.60mm
Overall Height	1016.0mm
Turning Radius	660.40mm
Ground Clearance	76.200mm
Total Weight	142.88kg - 174.18kg
Max Speed	8.046720kph
Battery Size (2 req.)	34Ah / U1 or 55Ah / 22NF
Range Per Charge	up to 19.312km / 30.578km
Motor	2 - Pole / 240W
Controller	P&G VR2 - A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	152.40mm x 50.800mm
Drive Wheels	355.60mm x 76.200mm
Rear Wheels	152.40mm x 50.800mm
Standard Colors	Red, White, Blue & Gray

Metric Measurements

Product #	P335
Weight Capacity	300lb
Drivetrain	FWD
Seat Width	14in - 20in
Seat Depth	14in - 20in
Backrest Height	18in - 22in
Seat-to-Floor Height	17.5in - 19.5in
Seal Elevation (Opt.)	12in
Overall Length	39.6in
Overall Width	22.5in

Overall Height	41.5in
Turning Radius	26in without footplate
Ground Clearance	3in
Total Weight	350lb - 385lb
Max Speed	5.2mph
Battery Size (2 req.)	MK 60Ah, M34 Gel
Range Per Charge	up to 19mi
Motor	2 - Pole / 222W
Controller	R-NET 90A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	14in x 3in hybrid gel
Rear Wheels	9in hybrid gel
Standard Colors	Noir Black, Dove Grey, Sapphire Blue, Ruby Red & Camo Green

