

Ultra Specifications

Imperial Measurements

Product #	P325
Weight Capacity	136.08kg
Drivetrain	MWD
Seat Width	406.40mm - 558.80mm
Seat Depth	406.40mm - 558.80mm
Backrest Height	457.20mm - 711.20mm
Seat-to-Floor Height	Power Tilt or Pan: 419.10mm - 495.30mm; Captain: 508.00mm - 584.20mm
Seal Elevation (Opt.)	304.80mm
Overall Length	914.40mm
Overall Width	609.60mm
Overall Height	1016.0mm
Turning Radius	660.40mm
Ground Clearance	63.500mm
Suspension	Front and Rear
Total Weight	142.88kg - 174.18kg
Max Speed	8.046720kph
Battery Size (2 req.)	35Ah, U1; 55Ah, 22NF, 55Ah, 22NF Gel
Range Per Charge	up to 30.578km
Motor	2 - Pole / 240W
Controller	P&G VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	355.60mm x 76.200mm foam-filled
Rear Wheels	152.40mm x 50.800mm foam-filled
Standard Colors	White, Red, Blue and Charcoal Grey

Metric Measurements

Product #	P325
Weight Capacity	300lb
Drivetrain	MWD
Seat Width	16in - 22in
Seat Depth	16in - 22in
Backrest Height	18in - 28in
Seat-to-Floor Height	Power Tilt or Pan 16.5in - 19.5in; Captain: 20in - 23in
Seal Elevation (Opt.)	12in
Overall Length	36in
Overall Width	24in
Overall Height	40in

Turning Radius	26in
Ground Clearance	2.5in
Suspension	Front and Rear
Total Weight	315lb - 384lb
Max Speed	5mph
Battery Size (2 req.)	35Ah, U1; 55Ah, 22NF, 55Ah, 22NF Gel
Range Per Charge	up to 19mi
Motor	2 - Pole / 240W
Controller	P&G VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	14in x 3in foam-filled
Rear Wheels	6in x 2in foam-filled
Standard Colors	White, Red, Blue and Charcoal Grey



