

## Velocity Specifications

### Imperial Measurements

<b>Product #</b>	<b>P325</b>
<b>Weight Capacity</b>	136.08kg
<b>Drivetrain</b>	MWD
<b>Seat Width</b>	406.40mm - 558.80mm
<b>Seat Depth</b>	406.40mm - 558.80mm
<b>Backrest Height</b>	584.20mm
<b>Seat-to-Floor Height</b>	419.10mm - 495.30mm
<b>Seal Elevation (Opt.)</b>	304.80mm
<b>Overall Length</b>	914.40mm
<b>Overall Width</b>	609.60mm
<b>Overall Height</b>	1016.0mm
<b>Turning Radius</b>	660.40mm
<b>Ground Clearance</b>	63.500mm
<b>Total Weight</b>	145.15kg - 183.25kg
<b>Max Speed</b>	9.656064kph
<b>Battery Size (2 req.)</b>	55Ah, 22NF Gel, 55Ah, 22NF AGM
<b>Range Per Charge</b>	up to 28.968km
<b>Motor</b>	2 - Pole / 275W
<b>Controller</b>	P&G VR2 or R-NET
<b>Charger</b>	6 Amp
<b>Gradient</b>	6°
<b>Front Wheels</b>	152.40mm x 50.800mm - solid
<b>Drive Wheels</b>	355.60mm x 76.200mm foam-filled
<b>Rear Wheels</b>	152.40mm x 50.800mm - solid
<b>Standard Colors</b>	Green, Metallic Black, Pink, Metallic Royal Blue, Candy Apple Red or Platinum Silver

### Metric Measurements

<b>Product #</b>	<b>P325</b>
<b>Weight Capacity</b>	300lb
<b>Drivetrain</b>	MWD
<b>Seat Width</b>	16in - 22in
<b>Seat Depth</b>	16in - 22in
<b>Backrest Height</b>	23in
<b>Seat-to-Floor Height</b>	16.5in - 19.5in
<b>Seal Elevation (Opt.)</b>	12in
<b>Overall Length</b>	36in
<b>Overall Width</b>	24in
<b>Overall Height</b>	40in
<b>Turning Radius</b>	26in
<b>Ground Clearance</b>	2.5in

<b>Total Weight</b>	320lb - 404lb
<b>Max Speed</b>	6mph
<b>Battery Size (2 req.)</b>	55Ah, 22NF Gel, 55Ah, 22NF AGM
<b>Range Per Charge</b>	up to 18mi
<b>Motor</b>	2 - Pole / 275W
<b>Controller</b>	P&G VR2 or R-NET
<b>Charger</b>	6 Amp
<b>Gradient</b>	6°
<b>Front Wheels</b>	6in x 2in - solid
<b>Drive Wheels</b>	14in x 3in foam-filled
<b>Rear Wheels</b>	6in x 2in - solid
<b>Standard Colors</b>	Green, Metallic Black, Pink, Metallic Royal Blue, Candy Apple Red or Platinum Silver



