



**AXCEL** P335

FWD Powerbase



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [Sales@Meritshealth.uk](mailto:Sales@Meritshealth.uk)

# Features

- Independent hydraulic, high-performance suspension
- Large 9" casters to navigate obstacles and rough terrain
- Bright LED lights, included
- Tool-free, rear battery access
- Power recline adjustment 85° - 160°
- Power tilt adjustment 0° - 50°
- Adjustable 6" shear reduction recline
- Powered by MK batteries, rated up to 19 miles per charge

# Accessories

(Based on chair configuration)

- USB charger
- Oxygen tank holder
- Cup holder

# Options

- ASL alternative drive
- Therafin positioning products
- Power articulating footplate
- Push handles
- Proform headrest with EZ-Clamp bracket

# Specifications

Product#	P335
Weight Capacity	300 lbs.
Drive train	FWD
Seat Width	14-22"
Seat Depth	14-22"
Backrest Height	18-22"
Seat-to-Floor Height	17.5-19.5"
Seat Elevation (Optional)	12"
Overall Length	39.6"
Overall Width	24.5"
Overall Height	41.5"
Turning Radius	26" w/o footplate
Ground Clearance	3"
Total Weight	350-385 lbs.
Max Speed	6 mph
Battery Size (2 required)	MK 60Ah, M34 Gel
Range Per Charge	up to 19 mi.
Motor	2-Pole/222W
Controller	R-NET 120A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-Tipper
Drive Wheels	14"x 3" hybrid gel
Rear Wheels	9" hybrid gel
Color Options	<ul style="list-style-type: none"><li>■ Black Noir</li><li>■ Dove Gray</li><li>■ Sapphire Blue</li><li>■ Ruby Red</li><li>■ Camo Green</li></ul>



Features MK batteries



**merits** Enhancing life through mobility.

©Merits Health Products UK Ltd  
3 Haugh Road, Burntisland, Fife KY3 0BZ  
www.MeritsHealth.uk | +44 7443 876057 | Sales@MeritsHealth.uk



rev. 0125